Endo-Eze[®] AET[™]

Contra-Angle Handpiece

Instructions for Use

Endodontic Handpiece for Instrumentation of Root Canals
30° Oscillation, 4:1 Reduction, ISO E-Type Coupling, Push-Button Chuck
Please read these instructions before operating the handpiece

Description & Indications for Use

The Endo-Eze AET Contra-Angle Handpiece (REF/UP #1154 or REF/UP# 1155) is a reciprocating contra-angle handpiece with a unique back-and-forth milling action used for root canal preparation and endodontic therapy as part of the Endo-Eze AET system. The Endo-Eze AET system and all associated components comprise a unique technique for canal shaping and cleaning known as Anatomic Endodontic Technology (AET).

Symbols

	This symbol means that handpiece exterior can be rinsed and disinfected in thermal washer disinfectors. For interior cleaning, Ultradent approves thermal washer disinfectors manufactured by Miele. Following every thermal washer-disinfection procedure, ensure that the product is cleaned and prepared using handpiece cleaning products from Kavo, W&H, or similar.
	This symbol means that handpiece can be sterilized in an autoclave up to a maximum temperature of 135° C. Ensure that handpiece is fully cleaned and prepared before any sterilization procedure. Ensure that handpiece is taken out of sterilizer as soon as sterilization is complete.

Precautions & Warnings

- The Endo-Eze AET Contra-Angle Handpiece is designed to be used with Endo-Eze Shaping Files only.
- Must be cleaned, prepared and sterilized prior to initial use. See section entitled "Cleaning, Preparing and Sterilizing Endo-Eze AET Handpieces" for more information.
- May be used with any ISO E-type motor.
- When operating handpiece, always consider safety of patient.
- Follow Endo-Eze AET instructions for technique and procedure.
- Do not drop or hit the handpiece against a hard object.
- Do not disassemble or alter the handpiece.
- Before use, place handpiece onto E-type motor, run for a few seconds and wipe away any excess oil. Note: If handpiece does not function normally, cease operation and return to Ultradent for repair.

Cleaning, Preparing and Sterilization of Endo-Eze AET Handpieces

Precautions and Warnings for Cleaning, Preparing and Sterilization:

- To minimize risk of cross infection, ensure that protective gloves are worn throughout entire cleaning, disinfection and sterilization process.
- Ensure that handpiece is cleaned and lubricated before any sterilization procedure. See section entitled "Lubrication" for more information on handpiece lubrication.
- Sterilize handpiece after every patient.
- Sterilize only with distilled water.
- Ensure that handpiece is removed from sterilizer as soon as sterilization is complete.
- Do not use ultrasonic cleaning units to clean Endo-Eze contra-angle handpiece.
- It is recommended that handpiece be sterilized immediately after use.
- The fittings on handpiece should be cleaned and cared for at least once per week using KaVo spray cleaner or similar agent.

Step by Step Instructions for Cleaning, Preparing and Sterilization

1. Remove contra-angle handpiece from E-type motor.
2. If necessary, detach file from its chuck/fitting.
3. Clean and disinfect the exterior of the handpiece and head using alcohol-soaked gauze, cotton roll or cellulose paper.
4. Clean and disinfect the interior of the handpiece using thermal washer-disinfector.
5. Lubricate handpiece with KaVo spray cleaner or similar agent.
6. Sterilize handpiece in sterilization bags in a steam sterilizer (autoclave) up to a maximum temperature of 135° C for 15 minutes.
Ensure that handpiece is taken out of sterilizer as soon as sterilization is complete; otherwise, the instrument may malfunction.
7. Before re-using, place Endo-Eze AET contra-angle handpiece onto an E-type motor, run for a few seconds, and wipe away any excess oil and lubricant.

Lubrication

Precautions and Warnings for Lubrication

- Lubricate before every sterilization procedure or twice daily after 30 minutes of continuous use.
- Before commencing lubrication, remove handpiece from its tubing connection and the Endo-Eze shaping file from the chuck.
- Shake the lubricant can 3-4 times to mix lubricant and propellant.
- Hold handpiece firmly to prevent it from slipping off spray cap during lubrication.
- Keep lubricants and lubricated handpiece away from open flames or ignition sources.
- Do not smoke when using this device or at any time during a lubrication procedure.

Step by Step Instructions for Lubrication

1. Screw E-type nozzle onto the spray cap over the arrow-head nozzle. Turn approximately 10 times.
2. Insert spray cap onto spray can.
3. Shake spray can heavily several times.
4. Insert the nozzle into the back end of the contra-angle handpiece.
5. Keep the spray can upright.
6. Spray for approximately 1-2 seconds.

Ultradent Recommendations for Cleaning, Disinfecting and Sterilizing Handpiece

External cleaning	- manual: - machine:	- brush the instrument under running water - clean the exterior of the instrument using a 60% alcohol solution - Miele Dental Washer Disinfector G 7781 or compatible device
Internal cleaning	- machine:	- Miele Dental Washer Disinfector G 7781 or compatible device
External disinfection	- manual: - machine:	- KaVo ELASTOclean 2, S&M MicroCid or similar surface disinfectant - Miele Dental Washer Disinfector G 7781 or compatible device
Internal disinfection	- machine:	- Miele Dental Washer Disinfector G 7781 or compatible device
Lubrication	- manual:	- KaVo spray 2112 A, or similar lubricating agent
Sterilization	- in sterilization bag	- Steam sterilizers (autoclave class B or S) B class: 135° C for 5 minutes S class: 135° C for 12 minutes

Technical Data - Contra-angle Handpiece	REF/UP 1154 or REF/UP 1155
Transmission Ratio	4:1
Color Coding	Green ring
Motor coupling according to standard	ISO 3964
Recommended root canal files	Endo-Eze shaping files with shank of 2.35 mm
Minimum chucking length	Until fully located
Maximum rated reamer/file speed (min-1)	6,000
Spray water volume (mL/min) according to ISO 7785-2	REF/UP 1155: >50 REF/UP 1154: No water spray
Recommended water pressure	28-43 psi (2-3 bar)
Recommended chip air pressure psi (bar) must be higher than water pressure	28-56 psi (2-4 bar)

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In case of service or repair, contact Ultradent Products, Inc.
This device is restricted to sale by or on the order of a dentist.
Keep out of reach of children.

For immediate reorder and/or complete descriptions of Ultradent's product line, refer
to Ultradent's catalog or call Toll Free **1-800-552-5512**.
Outside U.S. call (801) 572-4200
or visit our website at www.ultradent.com.

© Copyright 2008 Ultradent Products, Inc. All Rights Reserved.
Manufactured by NSK on behalf of Ultradent Products Inc.

Distributed by Ultradent Products Inc.
505 West 10200 South
South Jordan, Utah 84095, USA
#12706.1 030508

 **ULTRADENT**
PRODUCTS, INC.

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Fabricado por NSK en representación de Ultradent Products Inc.
Para mantenimiento o reparación, comuníquese con Ultradent Products, Inc.