**Endo-Eze®**  
**AET™**  
E-Type Air Motor

**Description & Indications for Use**

The Endo-Eze AET Air Motor (REF/UP #1178) is a standard E-type motor. It is intended for use with the Endo-Eze AET endodontic system. The Endo-Eze AET system and all associated components comprise a unique technique for canal shaping and cleaning known as Anatomic Endodontic Technology (AET).

**The Endo-Eze AET E-Type Air Motor contains the following.**

1. Coupling, handpiece side
2. O-Rings – (3) three O-rings, spray connection integrated
3. Mechanical direction. Clockwise and counter-clockwise rotation direction selected with speed regulation ring.
  - a. F = Clockwise
  - b. R = Counter-clockwise
4. Supply hose connection threads
5. Rubber gasket

**Symbols**

This symbol means that air motor exterior can be rinsed and disinfected in thermal washer disinfectors. For interior cleaning, Ultradent approves thermal washer disinfectors manufactured by Miele. Following every thermal washer-disinfection procedure, ensure that the product is cleaned and prepared using air motor cleaning products from Kavo, W&H, or similar.

This symbol means that air motor can be sterilized in an autoclave up to a maximum temperature of 135° C. Ensure that the air motor is fully cleaned and prepared before any sterilization procedure. Ensure that air motor is taken out of sterilizer as soon as sterilization is complete.

**Precautions & Warnings**

1. The Endo-Eze AET Air Motor is designed to be used with the Endo-Eze AET contra-angle, but it can be used with any standard E-Type attachment.
2. When removing the air motor from the air line, make sure that the speed regulating ring is set to position F or R to prevent it coming off due to any buildup of pressure. Do not remove the air motor during operation.
3. Air motor must be cleaned, prepared and sterilized prior to initial use. See section entitled "Cleaning, Preparing and Sterilizing Endo-Eze AET Air Motor" for more information.
4. When operating air motor, always consider safety of patient.
5. Follow Endo-Eze AET instructions for technique and procedure.
6. Do not drop or hit the air motor against a hard object.
7. Do not disassemble or alter the air motor.
8. Before use, place the attachment you intend to use onto E-type motor, run for a few seconds and wipe away any excess oil. Note: If air motor does not function normally, cease operation and return to Ultradent for repair.

**I. Cleaning, Preparing, Lubricating and Sterilizing Endo-Eze AET Air-Motor**

*Precautions and Warnings for Cleaning, Preparing and Sterilizing*

- To minimize risk of cross infection, ensure that protective gloves are worn throughout entire cleaning, disinfection and sterilization process.
- Ensure that the air motor is cleaned and lubricated before any sterilization procedure. See section below entitled "Lubrication of Endo-Eze AET Air Motor" for more information on lubrication of the air motor.
- Sterilize air motor after every patient.
- Sterilize only with distilled water.
- Ensure that air motor is removed from sterilizer as soon as sterilization is complete.

*Step by Step Instructions for Cleaning, Preparing and Sterilization*

1. Remove Endo-Eze AET E-type air motor from the air line. When removing the air motor from the air line, make sure that the speed regulating ring is set to position F or R to prevent it coming off due to any buildup of pressure. Do not remove the air motor during operation.
2. Clean and disinfect the exterior of the air motor using alcohol-soaked gauze, cotton roll or cellulose paper.
3. Clean and disinfect the interior of the air motor using thermal washer-disinfectant.

- Lubricate air motor with KaVo spray cleaner or similar agent.
- Sterilize air motor in sterilization bags in a steam sterilizer (autoclave) up to a maximum temperature of 135° C for 15 minutes. Ensure that air motor is taken out of sterilizer as soon as sterilization is complete; otherwise, the device may malfunction.
- Before re-using air motor, place Endo-Eze AET contra-angle handpiece (or equivalent) onto the air motor, run for a few seconds, and wipe away any excess oil and lubricant.

## Ultraviolet Recommendations for Cleaning, Disinfecting and Sterilizing Air Motor

External cleaning	- manual:  - machine:	- brush the instrument under running water - clean the exterior of the instrument using a 60% alcohol solution  - Miele Dental Washer Disinfector G 7781 or compatible device
Internal cleaning	- machine:	- Miele Dental Washer Disinfector G 7781 or compatible device
External disinfection	- manual:  - machine:	- KaVo ELASTOclean 2, S&M MicroCid or similar surface disinfectant  - Miele Dental Washer Disinfector G 7781 or compatible device
Internal disinfection	- machine:	- Miele Dental Washer Disinfector G 7781 or compatible device
Lubrication	- manual:	- KaVospray 2112 A, or similar lubricating agent
Sterilization	- in sterilization bag	- Steam sterilizers (autoclave class B or S) B class: 135° C for 5 minutes S class: 135° C for 12 minutes

### Lubrication of Endo-Eze AET E-Type Air Motor

#### Step by Step Instructions for Lubrication of Endo-Eze AET E-Type Air Motor

- Insert spray cap (with nozzle) onto spray can. Remove blue insert in nozzle and replace with black cone-shaped insert.
- Shake spray can several times.
- Set the speed ring of the air motor to maximum position (position F or R).
- Fit the drive air tube (the smallest air intake) of motor onto the spray cap.
- Keep spray can upright.
- Spray for approximately one second.

### II. Using the Endo-Eze AET Air Motor and Contra-Angle

#### Operating Pressure

- The operating pressure is measured as flow pressure.
- Operating pressure possible with 32 psi (2.2 bar) to 43 psi (3.0 bar).
- If only the air motor is used on the supply hose, set the pressure to 43 psi (3.0 bar).
- If the air motor is used alternately with a turbine on the same supply hose, it is essential to ensure that the pressure set for the two drives is within the prescribed range.

#### Operating Instructions

- Connect air motor to air supply hose. Align the motor back-end tubes with holes in the supply hose connector.
- Push the motor onto the supply hose connector.
- Turn the supply hose nut and securely fasten it. If loose, air and/or water may leak, and air pressure to the motor will become insufficient.

CAUTION: When motor speed feels low, check this connection first.

4. Firmly seat contra-angle onto air motor until it clicks into place.
5. Test run the air motor.
  - a. Start the air motor. Screw speed regulation ring to "F" clockwise rotation or "R" counter-clockwise rotation.
  - b. Vary the air motor from minimum to maximum within 30 seconds.
6. Test run the contra angle.
  - a. Hold contra-angle with head downward and let run approximately 30 seconds (without root canal file attached.)
7. To remove the contra-angle from the air motor, press the push-button on the air motor.
8. To remove the air motor, unscrew the supply hose nut and remove the air motor by pulling it away from the supply hose.

WARNING: When removing the air motor, make sure that the speed regulating ring is set to position F or R to prevent it coming off due to any buildup of pressure. Do not remove the air motor during operation.

#### *Changing the O-Rings*

- O-Rings can be exchanged without a special tool.
- Do not use a sharp tool.
- Firmly squeeze together the worn O-ring between the thumb and forefinger so as to form a loop and then withdraw the O-ring.

WARNING: Exchange damaged O-rings immediately!

Technical Data - Endo-Eze AET Air Motor	REF/UP 1178
Rated Speed (RPM) at 32 psi (2.2 bar)	5,000 to 20,000
Rated Speed (RPM) at 43 psi (3.0 bar)	5,000 to 25,000
Stationary torque (Ncm)	4
Maximum power (W)	20
Air consumption (Nl/min) at 32 psi (2.2 bar)	42
Air consumption (Nl/min) at 43 psi (3.0 bar)	50
Spray water volume (mL/min)	> 100
Operating pressure psi (bar)	32 psi to 43 psi (2.2 to 3.3 bar)
Exhaust air pressure psi (bar)	7.15 psi (<0.5 bar)
Water pressure psi (bar)	21.8 psi to 36.3 psi (1.5 to 2.5 bar)
Spray air pressure psi (bar) must be higher than the water pressure	28 psi to 43 psi (2 to 3 bar)

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In case of service or repair, contact Ultradent Products, Inc.  
This device is restricted to sale by or on the order of a dentist.  
Keep out of reach of children.

For immediate reorder and/or complete descriptions of Ultradent's product line, refer  
to Ultradent's catalog or call Toll Free **1-800-552-5512**.  
Outside U.S. call (801) 572-4200  
or visit our website at [www.ultradent.com](http://www.ultradent.com).

© Copyright 2008 Ultradent Products, Inc. All Rights Reserved.  
Manufactured by NSK on behalf of Ultradent Products Inc.

---

Distributed by Ultradent Products Inc.  
505 West 10200 South  
South Jordan, Utah 84095, USA  
#12705.1 030508

 **ULTRADENT**  
PRODUCTS, INC.

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Fabricado por NSK en representación de Ultradent Products Inc.
Para mantenimiento o reparación, comuníquese con Ultradent Products, Inc.